



Background

The RSV life-cycle within infected airway epithelial cells utilises positive sense (+ve sense) RNA replicative intermediates (RIs) including antigenome and mRNA transcripts for replication and viral protein synthesis respectively. hVIVO have now developed a strand-specific RT-qPCR assay to quantify specifically the active intra-cellular +ve sense RSV RI present in successfully infected epithelial cells. We demonstrate specific detection of <0.1 % +ve sense RI RSV in 99.9% virion genomic RSV mixtures with a lower limit of quantitation of RI RSV at 10 copies/assay well. Furthermore we have used the strand specific RI RSV RT-qPCR assay to generate area-under-the-curve (AUC) profiles using the cellular component of nasal washes collected from infected human volunteers and significantly different profiles of RI RSV AUC profile compared to total RSV AUC. Furthermore, the RI RSV AUC closely resembles infectious viral AUC generated by cell-based plaque assay data.

Introduction




Figure 1: RSV replicative cycle





Figure 2: Nasal wash samples from RSV inoculated volunteers were analysed by RT-qPCR (solid line, filled triangles), and plaque assay (broken line, empty circles) and plotted against collection time. (shaded towers equals symptoms score.)

RSV replicates through +ve sense stand replicative RNA intermediates (RI) in infected patients (Figure 1).

The area-under-the curve (AUC) data for RSV replication in RSV infected human volunteers is measured by total RT-qPCR detecting mainly extra-cellular genomic RNA (-ve strand). It is well-known that RSV AUC RT-qPCR curve has an exaggerated increase at late infection compared to plaque assay AUC data (Figure 2), This is caused by accumulation of non-inhibitable, replicative-defective particles. This late-stage AUC surge confounds antiviral efficacy determination of test compounds targeting RSV replication.

We hypothesise that measurement of intra-cellular RI RSV by RT-qPCR gives a fast, sensitive measure of anti-viral efficacy.

RI-RNA strand specific RT-qPCR assay in RSV infected humans



Nasal wash samples were collected from 4 RSV Memphis 37B inoculated healthy volunteers across study day 2 and day 12. Each plot is identified by the subject ID.

SUPERNATANT FRACTION

Total RT-qPCR (dotted line, □ symbol) and plaque assay data (solid line, ▼ symbols).

CELLULAR FRACTION

Total RT-qPCR (solid line, ▲ symbols) and RI specific RSV. (solid red line)

Both total RT-qPCR profiles (supernatant and cellular) show the characteristic high AUC late in infection (post day 8).

The RI-RSV AUC profile show a completely different profile with a lower maximum peak and rapid "cliff-face" fall-off late in infection. This profile is very similar to the plaque assay data.

Strand Specific reverse transcription (RT)-qPCR




Figure 3 Scheme for strand specific RT-qPCR

Mixture	RI-RSV copies	Genomic RSV copies	Ratio RI/gen
1	10 ¹	10 ⁴	1/1000
2	10 ²	10 ⁴	1/100
3	10 ³	10 ⁴	1/10
4	10 ⁴	10 ⁴	1/1
5	10 ⁵	10 ⁴	10/1
6	10 ⁶	10 ⁴	100/1
7	10 ⁷	10 ⁴	1000/1




Figure 4: Strand specific RT-qPCR was performed on in vitro RNA mixtures at different ratios (see inset table).

In vitro RNA copies of the RSV target region were synthesised in both +ve sense (RI) and -ve sense (genomic) orientations. Mixtures were made to generate 0.1% RI in genomic RNA background to 0.1% genomic in RI background (see mixture columns). Strand-specific RT-qPCR was performed on all the mixtures (see Figure 4).

Note the Ct difference between the 0.1% RI and genomic is Ct = 9, as expected for 3 log concentration change in template demonstrating 100% SPECIFICITY STRAND DETECTION.

The lower limit of quantitation for the RI-RSV RT-qPCR assay is 10 copies/rxn

CONCLUSION

1. TWO independent anti-viral efficacy measurements from ONE clinical sample is possible.

Cell pellet → RI-RSV specific RT-qPCR
Supernatant → total RT-qPCR (totRSV).

2. RI-RSV AUC profile is very similar to plaque assay.

Replicative/mRNA RSV has short t_{1/2} and rapid turnover whereas the genomic RSV in extracellular virion is stable. The totRSV AUC in late stage infection likely due to accumulation of replication-incompetent or neutralised with bound IgG. Either way, the intra-cellular RI-RSV AUC appears to reflect behaviour of of replication competent virions plaque-forming RSV.

3. RI-RSV AUC metric for virus replication inhibitors.

Unlike the plaque assay; no confounding cytotoxic IMP effects will be observed in RT-qPCR format. The RI-RSV AUC may be a more sensitive metric for anti-viral efficacy. The RI-RSV assay has a faster turnaround time compared to plaque assay